

Application Note - Digestion of various food matrices

Introduction:

This study evaluated the effectiveness of the **NOVAWAVE** microwave digestion system to digest different food groups like high fat products, proteins, vegetables and milk base products.

Sample Type:

CRM foods

- Sample weight: 0.5 g - 2 g
- 2 replicates

Supplies and Reagent:

- 1) **NOVAWAVE** model SA Microwave Digestion System*
- 2) Quartz 75 ml vessels, calibrated at 50 ml*
- 3) Teflon® Caps and Safety Pressure Release Caps, pre-set release pressure 30 bar (435 psi)*
- 4) Nitric Acid, **PlasmaPURE** HNO₃, 70%*
- 5) Analytical Balance 4 Decimal Places, Mettler-Toledo
- 6) ICP-MS PerkinElmer ELAN 6100
- 7) Spectroflame Modula FMD-07 ICP-OES, Spectro Analytical
- 8) Mini X-Flow Nebulizer*
- 9) 1.2mm Alumina Injector Torch*
- 10) Baffled Cyclonic Spray Chamber*
- 11) 10 ml Graduated Cylinder, Coming
- 12) **DigiTUBE**s*
- 13) **DigiFILTER**s*

*Manufactured by SCP SCIENCE.

Sample Preparation Procedure:

The samples were weighed on a 4 place analytical balance directly in the 75 quartz vessels. After adding 10ml of HNO₃, the samples were allowed to sit at room temperature for 5-10 minutes. The vessels were placed in the rack, capped and digested following the **NOVAWAVE** heating profile. After cooling to room temperature, the rack was vented carefully in a fumehood. The samples were normalized to 50 ml with deionized water. The oyster tissue and the toasted oat cereal were filtered prior to analysis.

NOVAWAVE Heating Program:

STAGE	RAMP TIME (Minutes)	PRESSURE (psi-Limit)	TEMPERATURE (°C)	HOLD TIME (Minutes)
1	10	435	190	10

*Note: This **NOVAWAVE** heating program is applied to the following CRM: peanut butter, baking chocolate, meat homogenate, protein powder, slurred spinach and infant formula. Only certified values are shown.

SCP SCIENCE

Providing Innovative Solutions to Analytical Chemists

NOVAWAVE

A new category of Automated Microwave Digestion

Application Note - Digestion of various food matrices

Results:

Peanut Butter (CRM 2387):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	435	75	411	18	106
Cu	324.7	4.7	0.1	4.93	0.15	96
Fe	259.9	16.4	1.1	16.4	0.8	100
K	766.5	5644	625	6070	200	93
Mg	279.5	1796	202	1680	70	106
Mn	257.6	17.1	1.9	16.0	0.6	106
Na	589.6	5393	279	4890	140	110
P	253.6	3156	53	3378	92	93
Zn	213.1	28.6	3.9	26.3	1.1	109

Baking Chocolate (CRM 2384):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	896	82	840	74	107
Fe	259.9	122	9	132	11	92

Meat Homogenate* (CRM 1546):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	330	3	323	28	102
Fe	259.9	12.2	0.8	11.4	1.0	107
Na	589.6	10783	465	9990	716	108

*Definition: a mixture of pork and chicken products blended together in a commercial process

Protein Powder (CRM 3244):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	12541	743	13280	900	94
Cu	324.7	11.7	1.2	10.2	1.0	114
K	766.5	17416	792	16000	1800	109
Mg	279.5	3377	188	3100	120	108
Mn	257.6	31.7	2.7	30.0	1.4	106
Na	589.6	896	65	910	100	98
P	253.6	12937	572	12200	880	93
Zn	213.1	120.7	1.4	126.4	7.7	95

NOVAWAVE

A new category of Automated Microwave Digestion

Application Note - Digestion of various food matrices

Slurried Spinach (CRM 2385):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	597	44	624	40	96
Fe	259.9	17.3	0.9	17.1	1.9	101
K	766.5	3377	145	3650	250	103
Mg	279.5	371	25	368	30	101
Mn	257.6	3.66	0.11	3.81	0.10	96
P	253.6	333.5	16.6	323.7	6.6	103
Zn	213.1	8.27	0.43	8.37	0.37	99

Infant Formula (CRM 1849):

Element Symbol	λ (nm)	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Ca	422.6	4628	261	4900	130	94
Cu	324.7	19.03	0.57	20.29	0.43	94
Fe	259.9	177.2	12.2	177.1	3.3	100
K	766.5	8762	342	8860	130	98
Mg	279.5	1667	38	1578	69	105
Mn	257.6	49.75	0.21	51.00	0.53	97
Na	589.6	4403	384	4150	140	106
P	253.6	3993	146	3782	36	105
Zn	213.1	162.3	4.4	152.3	5.1	106

Application Note - Digestion of various food matrices

Oyster Tissue (CRM 1566b)

NOVAWAVE Heating Program:

STAGE	RAMP TIME (Minutes)	PRESSURE (psi-Limit)	TEMPERATURE (°C)	HOLD TIME (Minutes)
1	10	435	210	20

Element Symbol	M/Z amu	NOVAWAVE Conc (ppm)	NOVAWAVE SD (ppm)	Cert. Value Conc (ppm)	Cert. Value SD (ppm)	Recovery (%)
Na	23	3286	102	3297	53	99
Mg	24	1103	73	1085	23	102
Al	27	185.8	14.5	197.2	6.0	94
K	39	6489	464	6520	90	99
Ca	44	874	64	838	20	104
V	51	0.516	0.051	0.577	0.023	89
Mn	55	19.6	0.9	18.5	0.2	106
Fe	57	221.1	12.4	205.8	6.8	107
Co	59	0.401	0.02	0.371	0.009	109
Ni	60	0.93	0.07	1.04	0.09	90
Cu	63	76.6	11.2	71.6	1.6	108
Zn	66	1529	172	1424	46	107
Se	78	2.21	0.24	2.06	0.15	107
Rb	85	3.616	0.706	3.262	0.145	110
Ag	107	0.616	0.02	0.666	0.009	92
Cd	111	2.42	0.30	2.48	0.08	98
Pb	206	0.331	0.026	0.308	0.009	108

Application Note - Digestion of various food matrices

Salad Dressing

NOVAWAVE Heating Program:

STAGE	RAMP TIME (Minutes)	PRESSURE (psi-Limit)	TEMPERATURE (°C)	HOLD TIME (Minutes)
1	10	435	200	10

*Note: This is not a NIST CRM

Element Symbol	M/Z amu	NOVAWAVE Conc (ppm)	Known Value Conc (ppm)	Recovery (%)
Ca	44	544	548	99
Cu	63	0.44	0.43	102
Fe	57	1.99	1.89	106
K	39	1722	1590	108
Mg	24	88.8	83.3	107
Mn	55	0.26	0.28	93
Na	23	7273	7540	96
P	31	489	454	108
Zn	66	2.4	2.2	109

Toasted Cereal Oats

A) Sample Size: 0.5 g

NOVAWAVE Heating Program:

STAGE	RAMP TIME (Minutes)	PRESSURE (psi-Limit)	TEMPERATURE (°C)	HOLD TIME (Minutes)
1	10	435	200	10

*Note: This is not a NIST CRM

Element Symbol	λ nm	NOVAWAVE Conc (ppm)	Known Value Conc (ppm)	Recovery (%)
Ca	422.6	4564	4290	106
Cu	324.7	2.60	2.73	95
Fe	259.9	249	241	103
K	766.5	2414	2440	99
Mg	279.5	911	860	106
Mn	257.6	21.2	22.4	95
Na	589.6	9489	9670	98
P	253.6	3586	3390	106
Zn	213.1	218	203	107

Application Note - Digestion of various food matrices

B) Sample Size: 2 g

NOVAWAVE Heating Program:

STAGE	RAMP TIME (Minutes)	PRESSURE (psi-Limit)	TEMPERATURE (°C)	HOLD TIME (Minutes)
1	20	435	210	25

Element Symbol	λ nm	NOVAWAVE Conc (ppm)	Known Value Conc (ppm)	Recovery (%)
Ca	422.6	4039	4290	94
Cu	324.7	2.88	2.73	106
Fe	259.9	250	241	104
K	766.5	2550	2440	104
Mg	279.5	815	860	95
Mn	257.6	23.0	22.4	103
Na	589.6	9110	9670	94
P	253.6	3370	3390	99
Zn	213.1	223	203	110

Discussion / Conclusion:

After the digestion and subsequent dilution, most solutions are clear and no sediments are noticeable. The oyster tissue and the toasted oat cereal showed some sedimentation but did not affect the recovery of the results. The recoveries are within 10% of the NIST certified values and/or within the standard deviations of the published NIST value. Hence, the **NOVAWAVE** is appropriate for digesting a variety of food matrices.

The **NOVAWAVE** has a dynamic temperature range allows for many types of food to be completely digested with just nitric acid. The ICP-OES (MS) analysis is based upon the best performance across the wavelength (atomic mass) range required to analyze all of the elements at once. The calibration curves are typically over 3 orders of magnitude with each CRM analyzed on one calibration curve.

A comparison was made using toasted cereal oats between 0.5 g and 2 g sample size to evaluate the effectiveness of a complete digestion with a sample size 4 times larger than typical employed. The recoveries of both sample sizes were comparable.

By:
René-Viet Nguyen, Ph.D, P. Chem.
Application Chemist

SCP SCIENCE

Providing Innovative Solutions to Analytical Chemists

www.scpscience.com
sales@scpscience.com

CANADA / USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

INTERNATIONAL

Tel.: +1 (514) 457-0701
Fax: +1 (514) 457-4499

EUROPE

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

CHINA

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302

BRAZIL

Tel.: +55 19 8194-8488