



INORGANIC

UNDERGRADUATE EXPERIMENT

Measuring the Isotopic Ratio of $^{10}\text{B}/^{11}\text{B}$ by 60 MHz ^1H NMR Spectroscopy




Via Guido Rossa 14
23875 Osnago (LC) - Italy
Tel. +39 039 9280061
Fax. +39 039 9289636
info@quantanalitica.com
www.quantanalitica.com


www.nanalysis.com

INTRODUCTION

NMR spectroscopy is used widely across all branches of chemistry due to its powerful structure elucidation capabilities and the inherently quantitative nature of this technique. The technique is typically introduced in the organic chemistry curriculum with discussion primarily focused on nuclei with a nuclear spin quantum number, I , of $\frac{1}{2}$, i.e. ^1H and ^{13}C . Typically, quadrupolar nuclei are specialized and require more advanced probes to detect. Therefore, they are often reserved for more advanced inorganic courses.



In this experiment, we employ a simple indirect method for an undergraduate student to accurately measure the distribution of boron isotopes by detecting the hydrogen atoms directly bonded to the boron atoms.^[1] While the boron isotopes are NMR active and could be observed at different resonant frequencies (see table 1), it is not possible to record spectra for both isotopes in a single 1D NMR experiment. While both ^{10}B and ^{11}B NMR experiments have been performed and studied separately,^[2] it is only through coupling to a common nucleus that both nuclei can be observed indirectly in a single spectrum. Therefore, the ^1H NMR spectrum of sodium borohydride is recorded and the effects of the boron nuclei on the resulting spectrum will be used as an indirect observation probe to calculate the isotopic ratio of $^{10}\text{B}/^{11}\text{B}$.



BACKGROUND

Crucial to the success of this experiment is the fact that ^{10}B and ^{11}B are both NMR active and have different nuclear spin quantum numbers (table 1). Therefore, each boron isotope will split the ^1H nuclei in BH_4^- differently.

The expected splitting pattern can be calculated using the

$$2n+1 \text{ rule}$$

where I = spin quantum number
 n = number of nuclei

Students are likely familiar with the $n+1$ rule used to predict the splitting pattern in a ^1H NMR spectrum, which is derived from the $2n+1$ rule. Since I is $\frac{1}{2}$ for the ^1H nucleus, the equation simplifies into $n+1$.

Table 1. Properties of Selected NMR-Active Nuclei^[3]

Isotope	I	Gyromagnetic Ratio (γ) ^a	% Natural Abundance	Larmor Frequency ^b
^1H	$\frac{1}{2}$	26.7522	99.99	400.00
^{13}C	$\frac{1}{2}$	6.7282	1.07	100.58
^{10}B	3	2.8747	19.90	42.98
^{11}B	$\frac{3}{2}$	8.5847	80.10	128.34


^aGyromagnetic ratios are in units of $10^7 \text{ rad s}^{-1} \text{ T}^{-1}$.

^bLarmor frequencies in MHz at a magnetic field of 9.4 T.

For this particular example, n is 1 for both cases as there is only one boron nuclei in the borohydride molecule. Therefore, the ^{10}B isotope, with $I = 3$, will split the coupled ^1H nucleus into a heptet ($2 \times 3 \times 1 + 1$) while ^{11}B , which has $I = \frac{3}{2}$, will result in a quartet ($2 \times \frac{3}{2} \times 1 + 1$). Due to the distinctly different coupling constants between a proton and the two boron isotopes, all 11 peaks are well resolved and integrating the individual signals are straightforward. Comparing the sum of the areas of the heptet (^{10}B) with the combined areas of the quartet (^{11}B) will give the isotopic distribution of the two boron isotopes.

PROCEDURE

Note: The cap should be vented to prevent pressure build-up due to the slow production of hydrogen gas as the NaBH_4 slowly reacts with D_2O .



RESULTS



DISCUSSION

Given the $2I_n + 1$ rule:

$^{11}\text{B}: (2 \times (3/2) \times 1) + 1 = 4$

$^{10}\text{B}: (2 \times 3 \times 1) + 1 = 7$

Because the two signals overlap, instead of summing the multiplets traditionally, we integrate each peak individually and sum them.

Total integration = 3.992 + 0.897 = 4.889
 $\%^{11}\text{B} = (3.992/4.889) \times 100 = 81.65\%$
 $\%^{10}\text{B} = (0.897/4.889) \times 100 = 18.34\%$

CONCLUSIONS

REFERENCES

- ## DATA ACCESSIBILITY

- 1) Chemosensors and ^{19}F NMR Spectroscopy
- 2) Isomerization of Mo complexes via ^{31}P NMR Spectroscopy
- 3) Aldol Condensation

www.nanalysis.com/sample-experiments

 **www.nanalysis.com**
 **1.855.NMReady**
 **@nanalysis**

