



Range spettrale da 220 a 3800 nm



CYLINDRICAL CELLS



Catalog Number	Description	Path Length	Volume (ml)
C-4010	Cylindrical cell w/PTFE stopper (1)	10mm	3
C-4020	Cylindrical cell w/PTFE stopper (1)	20mm	6
C-4040	Cylindrical cell w/PTFE stopper (1)	40mm	12
C-4050	Cylindrical cell w/PTFE stopper (2)	50mm	15
C-4100	Cylindrical cell w/PTFE stopper (2)	100mm	30
C-4010-T	Cylindrical cell w/threaded top & screw cap (1)	10mm	3
C-4020-T	Cylindrical cell w/threaded top & screw cap (1)	20mm	6
C-4040-T	Cylindrical cell w/threaded top & screw cap (1)	40mm	12
C-4050-T	Cylindrical cell w/threaded top & screw cap (2)	50mm	15
C-4100-T	Cylindrical cell w/threaded top & screw cap (2)	100mm	30
C-4010-G	Cylindrical cell w/glass tube attached (1)	10mm	3
C-4020-G	Cylindrical cell w/glass tube attached (1)	20mm	6
C-4040-G	Cylindrical cell w/glass tube attached (1)	40mm	12
C-4050-G	Cylindrical cell w/glass tube attached (1)	50mm	15
C-4080-G	Cylindrical gas cell w/quartz to borosilicate graded seal tube attached (1)	80mm	24
C-4100-G	Cylindrical cell w/glass tube attached (1)	100mm	30

LARGE DIAMETER CYLINDRICAL CELLS

Catalog Number	Description	Path Length	Volume (ml)
LDC-4010-T	Cylindrical cell - 50 mm outside diameter, w/threaded top & screw cap (1)	10mm	17
LDC-4020-T	Cylindrical cell - 50 mm outside diameter, w/threaded top & screw cap (1)	20mm	35
LDC-4040-T	Cylindrical cell - 50 mm outside diameter, w/threaded top & screw cap (1)	40mm	70
LDC-4050-T	Cylindrical cell - 50 mm outside diameter, w/threaded top & screw cap (2)	50mm	86
LDC-4100-T	Cylindrical cell - 50 mm outside diameter, w/threaded top & screw cap (2)	100mm	172

SMALL DIAMETER CYLINDRICAL CELLS

Catalog Number	Description	Path Length	Volume (ml)
CRD-4050-T	Cylindrical cell - 10 mm outside diameter, w/threaded top & screw cap (1)	50mm	7
CRD-4100-T	Cylindrical cell - 10 mm outside diameter, w/threaded top & screw cap (1)	100mm	14



MICRO CYLINDRICAL CELLS

Catalog Number	Description	Path Length	Volume (ml)
MC-4005	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	0.05mm	1.9
MC-401	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	0.1mm	2.0
MC-402	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	0.2mm	2.1
MC-405	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	0.5mm	2.2
MC-410	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	1mm	2.3
MC-420	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	2mm	2.5
MC-450	Micro cylindrical cell - 22 mm outside diameter, w/PTFE stopper (1)	5mm	3.0

