



Range spettrale da 220 nm a 3800 nm



FLUOROMETER CELLS



Catalog Number	Description	Path Length	Volume (ml)
RF-4005	Open top w/PTFE lid, 5 polished windows	5x10	1.7
RF-4010	Open top w/PTFE lid, 5 polished windows	10x10	3.5
RF-4005-T	Threaded top w/PTFE screw cap, 5 polished windows	5x10	1.7
RF-4010-T	Threaded top w/PTFE screw cap, 5 polished windows	10x10	3.5
RF-4010-S	Stopper top w/PTFE stopper, 5 polished windows	10x10	3.5
RF-4005-G	Cell w/quartz to borosilicate graded seal tube attached, 5 windows	5x10	3.5
RF-4010-G	Cell w/quartz to borosilicate graded seal tube attached, 5 windows	10x10	3.5
RF-4010-P	Quartz tube attached, 5 polished windows	10x10	3.5
RF-4005-F	Flow cell with two tubes attached, open both ends, 4 polished windows	5x10	1.7
RF-4010-F	Flow cell with two tubes attached, open both ends, 4 polished windows	10x10	3.5
RF-4020	Open top w/PTFE lid, 5 polished windows	10x20	7.0
RF-4040	Open top w/PTFE lid, 5 polished windows	10x40	10.5

NOTE



RF-403-45

MINI FLUOROMETER CELLS




A-6305


Catalog Number	Description	Path Length	Volume (ml)
RF-403-45	Open top w/PTFE lid, 5 polished windows	3x3	0.3
RF-404-45	Open top w/PTFE lid, 5 polished windows	4x4	0.5
RF-405-45	Open top w/PTFE lid, 5 polished windows	5x5	0.8

CELL HOLDERS FOR MINI FLUOROMETER CELLS

A-6303	Cell Holder for 3x3x45mm (for RF-403)
A-6304	Cell Holder for 4x4x45mm (for RF-404)
A-6305 <i>(shown)</i>	Cell Holder for 5x5x45mm (for RF-405)



MIRO and SEMI-MICRO FLUOROMETER CELLS



Catalog Number	Description	Path Length	Volume (ml)
SMF-4010	Open top semi-micro w/PTFE lid, 5 polished windows	4x10	1.4
SMF-4010-T	Threaded top semi-micro w/PTFE screw cap, 5 polished windows	4x10	1.4
SMF-4010-G	Semi-micro cell w/quartz to borosilicate graded seal tube, 5 windows	4x10	1.4
SMF-4010-F	Semi-micro cell w/quartz to borosilicate graded seal tube, 5 windows	4x10	1.4
MLF-4010	Open top micro cell w/PTFE lid, 5 polished windows	2x10	0.6
MLF-4010-T	Threaded top micro cell w/PTFE screw cap, 5 polished windows	2x10	0.6
MLF-4010-G	Micro cell w/quartz to borosilicate graded seal tube attached, 5 windows	2x10	0.6
MLF-4010-F	Micro cell w/quartz to borosilicate graded seal tube attached, 5 windows	2x10	0.6