

Designed for Larger Volumes

Savillex pressure vessels are heavy-walled, PFA containers designed to accommodate larger volumes of material in high temperature and pressure applications. These containers can also be used to transfer large quantities of chemicals, or they can be utilized with our column components for a variety of other applications. Savillex offers three standard volumes and we can make larger volume units upon request.


In addition to our standard pressure vessel product offering, we also have the in-house capability to make custom configured pressure vessels to meet your specific volumes/lengths and fitting needs.

- Pressure vessels are 4.20" OD x 3.70" ID x customer defined length
- 75 PSI pressure rating up to 100° F
- Fittings are 1/4" pipe and 1/2" pipe sizes
- Welded fitting adaptors available in tubing flare, MNPT, FNPT

Please see our custom capabilities and contact us to discuss your specific needs.



End Cap Options



Custom Pressure Vessel Examples

