



Elemental Microanalysis

Now offering ICP-MS cones




OEM No.	EMAL No.	Description	Material	Suitable for the following Agilent Technologies® instruments						
				HP4500	7500a	7500c	7500cs	7500s	7500ce	7500cx
G1820-65050	V1000	Skimmer Cone	Nickel		x	x		x		
G3270-65024	V1001	Skimmer Cone	Nickel		x				x	x
G1820-65207	V1002	Skimmer Cone	Platinum		x					
G1833-65132	V1003	Skimmer Cone (Copper Base)	Platinum				x		x	x
G1820-65238	V1004	Sample Cone	Nickel	x	x	x	x	x	x	x
G1820-65239	V1005	Sample Cone (10mm)	Platinum	x	x	x	x	x	x	x
G1820-65265	V1006	Skimmer Cone	Nickel	x	x			x		
G1820-65237	V1007	Skimmer Cone (7.2mm)	Platinum	x	x			x		
G1833-65497	V1008	Skimmer Cone	Nickel			x	x			
G3270 -60106	V1009	Skimmer Cone (Nickel Base)	Platinum				x			

OEM No.	EMAL No.	Description	Material	Suitable for the following Agilent Technologies® instruments				
				7700x/e	7800	7900	8800	8900
G3280-67040	V1011	Sample Cone (Copper Base)	Nickel	x	x	x	x	x
G3280-67041	V1012	Skimmer Cone	Nickel	x	x		x	
G3280-67036	V1013	Sample Cone	Platinum	x	x	x	x	x
G8400-67200	V1014	Skimmer Cone	Nickel			x		x
G3280-67064	V1015	Skimmer Cone	Platinum	x		x	x	
G3280-67060	V1016	Skimmer Cone	Platinum	x	x		x	
G8400-67202	V1017	Skimmer Cone	Platinum			x		x
G3280-67066	V1018	Skimmer Cone	Nickel	x		x	x	
G3280-67063	V1019	Skimmer Cone (Nickel Base)	Platinum	x	x		x	
G8400-67201	V1020	X Skimmer Cone	Platinum			x		x
G3280-67065	V1021	Skimmer (Nickel Cone)	Platinum	x		x	x	
G3280-67061	V1022	'Plated' Sample Cone	Nickel	x	x	x	x	x
G3280-67056	V1023	Sampler Cone (18mm)	Platinum	x	x	x	x	x



Skimmer Cone (V1001)




Sample Cone (Copper Base) (V1011)

OEM No.	EMAL No.	Description	Material	Suitable for the following PerkinElmer® instruments		
				ELAN 250	ELAN 500	ELAN 5000
WE01-3087	V3000	Sample Cone	Nickel	x	x	x
N812-2000	V3001	Skimmer Cone	Nickel	x	x	x

OEM No.	EMAL No.	Description	Material	Suitable for the following PerkinElmer® instruments			
				ELAN 6000	ELAN 6100	ELAN 9000	ELAN DRC
WE02-1140	V3002	Sample Cone	Nickel	x	x	x	x
WE02-1137	V3003	Skimmer Cone	Nickel	x	x	x	x
WE02-7802	V3004	Sample Cone	Platinum	x	x	x	x
WE02-7803	V3005	Skimmer Cone	Platinum	x	x	x	x
WE01-7142	V3006	Cone Insertion Tool		x			

OEM No.	EMAL No.	Description	Material	Suitable for the following PerkinElmer® instruments				
				NexION 300	NexION 350	NexION 1000	NexION 2000	NexION 5000
W103-3995	V3007	Hyper Skimmer Cone	Aluminium	x	x	x	x	
W102-6356	V3008	Skimmer Cone	Nickel	x	x	x	x	x
W103-3612	V3009	Sample Cone	Nickel	x	x	x	x	x
W102-6907	V3010	Skimmer Cone	Platinum	x	x	x	x	x
W103-3614	V3011	Sample Cone	Platinum	x	x	x	x	x



Sample Cone (V3009)



Skimmer Cone (V3008)

‘Price is what you pay, value is what you get’.

As the leading independent supplier of elemental analysis consumables, **we can offer substantial savings in comparison to the prices charged by the original supplier of your instrument.** But we know that you need more from a supplier than simply low prices. You need to buy with the confidence that your analysis will not be compromised by poor quality consumables and the knowledge that your supplier has all the resources to provide the full range of support that you expect.

OEM No.	EMAL No.	Description	Material	Plasma Class	Suitable for the following NU Instruments Ltd equipment							
					Plasma		Plasma II		Plasma 3		AttoM	
					(Std)	(ES)	(Std)	(ES)	(Std)	(ES)	(Std)	(ES)
325-294	V2000	Skimmer Cone – HS1-7	Nickel	Dry			x	x	x	x	x	x
319-497	V2001	Wide Angle Skimmer Cone	Nickel	Wet	x	x	x	x	x	x		
319-595	V2002	Wide Angle Skimmer Cone	Nickel	Wet	x	x	x	x	x	x	x	x
319-540	V2003	Skimmer Cone – AB6	Nickel	Dry	x		x		x		x	
319-285	V2004	Wide Angle Sample Cone	Nickel	Wet	x	x	x	x	x	x	x	x
319-541	V2005	Sample Cone – AB9	Nickel	Dry	x	x	x	x	x	x	x	x
319-284	V2006	Narrow Angle Skimmer Cone	Nickel	Wet	x	x						x
319-656	V2007	Sample Cone (Screw-In)	Platinum	Wet				x		x		x
319-594	V2008	Wide Angle Skimmer Cone	Platinum	Wet	x	x	x	x	x	x		
319-593	V2009	Wide Angle Sample Cone	Platinum	Wet	x	x	x	x			x	x
319-599	V2010	Skimmer Cone	Nickel	Dry	x							



Skimmer Cone – AB6 (V2003)



Skimmer Cone – HS1-7 (V2000)

Our promise to you

We will always treat our customers with honesty and integrity. If we do make a mistake, we don't try to hide it – we recognise it and put it right.

We won't sell you short measure, and we won't mislead you in any way. Our business is built on long-term relationships, as our customers will testify – many have been with us for more than 30 years.

OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments			
				PQ2	VG ExCell	X-SERIES	iCapQ/Qnova
3201101	V4000	Micro-Skimmer Cone	Platinum	x	x	x	
3200860	V4001	Micro-Skimmer Cone	Nickel	x	x		
3004661	V4002	Sample Cone	Nickel	x	x	x	
3601289	V4003	Sample Cone	Platinum			x	x
3600812	V4004	Sample Cone	Nickel			x	x
3600811	V4005	Skimmer Cone	Nickel			x	
3601320	V4006	Skimmer Cone Adapter	Nickel			x	



Skimmer Cone (V4005)



'H' Skimmer Cone (V4008)

OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments		
				ELEMENT	ELEMENT 2	NEPTUNE
1044530	V4007	Sample Cone	Nickel	x	x	x
1067600	V4008	'H' Skimmer Cone	Nickel	x	x	x
1067500	V4009	Sample Cone	Platinum	x	x	x
1067501	V4010	Sample Cone (w/Boron)	Platinum	x	x	x
1047460	V4011	'H' Skimmer Cone	Platinum	x	x	x
1047461	V4012	'H' Skimmer Cone (w/Boron)	Platinum	x	x	x
1142160	V4013	'X' Skimmer Cone	Nickel	x	x	x
1142150	V4014	'X' Skimmer Cone (Boron free)	Platinum	x	x	x
1260630	V4015	Jet Sample Cone	Nickel		x	x

Guaranteed quality

Our company's quality system is certified to ISO 9001:2015 by LRQA, we have an ISO/IEC 17025:2017 UKAS accredited Testing Laboratory, and we are an ISO/IEC 17034:2016 UKAS accredited Reference Material Producer. Our confidence in our products allows us to offer you a total satisfaction guarantee: if our product does not perform at least as well as the brand of product you usually use, you can return it for a full refund.


OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments						
				VG ExCell	X-SERIES	iCapQ/ Qnova	ELEMENT	ELEMENT 2	NEPTUNE	Apix
3200456	V4016	Plasmatrace Skimmer Cone	Nickel	x						
3201941	V4017	Skimmer Cone (7.2mm insert)	Platinum	x						
3201585	V4018	Skimmer Cone (0.4 Hole – Eclipse)	Nickel	x						
1067510 TS067016	V4019	Skimmer Cone (N)	Platinum				x	x	x	
1044540 TS067015	V4020	Skimmer Cone (N)	Nickel				x	x	x	
1044850 FS067017B	V4021	Skimmer Fitting Tool	Resin				x	x	x	
1142151	V4022	'X' Skimmer Cone (with Boron)	Platinum				x	x	x	
1184050	V4023	Sample Cone	Aluminium				x	x	x	
1184060	V4024	Skimmer Cone	Aluminium				x	x	x	
3600813	V4025	Skimmer Cone Adapter (w/4 screws)	Copper		x					
4600347	V4026	Xi Interface Kit			x					
3601219	V4029	Torch Shield	Silver		x					
3600814	V4030	Screws for Skimmer Cone Adaptor (pk/4)	Stainless Steel		x					
3601245	V4031	Skimmer Cone Adaptor	Copper Nickel		x					
1341430	V4032	Skimmer Cone (Hot/ Cold Plasma)	Platinum			x				
3202521	V4033	Sample Cone, Insert (6.35mm)	Platinum	x						
1600335	V4034	Sampler	Platinum		X					



Micro-Skimmer Cone (V4001)




Sample Cone (V4007)



Sample Cone (V4004)

Elemental Microanalysis holds the following accreditations:



All trademarks acknowledged.
Use of OEM name or trademark does not imply
products are of OEM origin, unless specified.
Copyright © 2022 Elemental Microanalysis